MiaCara°

Nube Box-Bed









SKII:

C02-002-01-S/SM/M/L, C02-002-02-S/SM/M/L, C02-002-03-S/SM/M/L, C02-002-04-S/SM/M/L

Description:

The popular MiaCara box bed is now also available in a high quality yet robust Italian fabric in a structured look in 4 earthy colors.

The execution of our Box-Bed is typical MiaCara: purist design combined with high functionality and quality.

It is the inner values that count: for the filling we use high quality foam made in Germany.

The side panels are soft, yet sturdy enough to lay against and rest the head on. For maximum comfort, the mattress is filled with viscoelastic foam. This material adapts perfectly to the body shape of the dog and provides optimum support for the spine and relieves the joints. Visco foam is permanently elastic and offers a unique pressure relief for exceptional lying comfort.

The mattress is reversible in the same look. The cover has an inner zipper and can be easily removed. For maximum hygiene the cover and all foam parts are machine washable. Quality Made in Germany.

Sizes, Dimensions & RRP:

S: 70 x 55 x 22 cm outer dimension (L x W x H)

RRP: 229,- Euro

SM: 80x65x24 cm Außenmaß (L x W x H)

RRP: 269,- Euro

M: 91 x 71 x 26 cm outer dimension (L x W x H)

RRP: 299,- Euro

L: 120 x 93 x 30 cm outer dimension (L x W x H)

RRP: 399.- Euro

Available Colors:

- Salt&Pepper
- Tobacco
- Smoke
- Rosé











Salt&Pepper Tobacco

Smoke

Rosé

Material:

- Inner cushion: cotton blend fabric 60% Cotton/ 40% Polyester
- Foam parts: 100% Polyurethan
- · Cover: 61% Polyester/39% Cotton

Foam properties:

- · stable and durable
- · climatising
- hypoallergenic
- · permanently elastic

Fabric properties:

- · Light Fastness: 5
- Pilling: 4/5
- high abrasion resistance: 40.000 Martindale
- 2,5% shrinkage

Care Instructions:

Inner Cushion Sides: 60° C gentle wash cycle











The cushion filling is made of extremely soft, open cell foam. To obtain the excellent properties of this material in the long run, the cushion should be maintained and cleaned accordingly.

The best cleaning procedure for the inner cushion is to laundry it with a commercially available mild detergents at temperatures up to 60° C. It is not recommended to use highly chlorinated and chemically aggressive washing additives. A brownish discoloration of the foam can not be avoided, but has no effect on the quality of the foam.

After washing and rinsing please remove as much water as possible through the spin cycle. The inner cushion can be best dried by air on a rack. It can also be tumble dried at 90° C, but please avoid over-drying.

Please plump up the inner cushion regularly so that it retains its fluffiness and shape.

Mattress: 30° C gentle wash cycle











Wet wash with a commercially available mild detergent at temperatures up to 30° C. It is not recommended to use highly chlorinated and chemically aggressive washing additives as this can damage the memory foam when drying.

After washing and rinsing please remove as much water as possible through the spin cycle. The mattress can be best dried by air on a rack. It can also be tumble dried at 90° C.

Cover: 30° C gentle wash cycle











For simple and uncomplicated cleaning in between, just vacuum the surface of the cover with an upholstery brush.